



Platinum Polishing Kit

Product Specifications



The platinum polishing kit includes:

- Polishing Lathe (Variable Speed, 650mm length)
- 240V AC/50Hz or 120V AC/60Hz
- Crucible and Mould Polishing Jigs
- Polishing Paste (3 & 8 Diamond Paste)
- Polishing Oil & Pump Pack
- 1000 Grit Paper
- Soft & Hard Felt Pads

Part No. 951405 – 230VAC, 50Hz, VDE style plug

Part No. 951406 – 230VAC, 50Hz, SAA style plug

Part No. 951408 – 230VAC, 50Hz, UL style plug

CE and CSA/UL Approved – Weight 23.5Kg

XRF Chemicals

xrf.chemicals@xrfscientific.com

Tel: +61 8 9244 9600
Fax: +61 8 9244 9611

XRF Labware

xrf.labware@xrfscientific.com

Victoria
Tel: +61 3 9408 4811
Fax: +61 3 9408 6811

Western Australia
Tel: +61 8 9240 8929
Fax: +61 8 9240 8930



XRF Technology

xrf.technology@xrfscientific.com

Victoria
Tel: +61 3 9720 6339
Fax: +61 3 9720 6412

Western Australia
Tel: +61 8 9274 7222
Fax: +61 8 9274 7244

©2011 XRF Scientific Ltd. All rights reserved. All trademarks are the property of XRF Scientific Ltd and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.



©2011 XRF Scientific Ltd. All rights reserved. All trademarks are the property of XRF Scientific Ltd and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.