



AUTOMATIC WEIGHING: FLUX

XRWEIGH

LIMS



The xrWeigh takes the hard work out of flux weighing, removing a repetitive task from your laboratory while improving accuracy and traceability.

Simple Operation – “One touch” operation, just turn it on! Attach flux bottle, set your target and let the xrWeigh do the rest

Easy Clean – Simple flux bottle replacement, easy access for cleaning, low maintenance

Reliable – Based on the existing single system architecture, simple modular design, extensively tested

Key features

- In built stability checking system
- Adjustable load receptor with overload detection
- Internal calibration weight for improved accuracy and repeatability
- Force compensation technology weighing system
- Simple touch screen operation for selection of target weight and range

Technical specifications

- Capacity: 1 Vial or 1 crucible
- Readability: 0.1 mg
- Accuracy ≤ 1 mg (speed dependent)
- Linearity: $\pm .2$ mg
- Speed > 90 vials per hour
- Dimensions: 420mm(W) x 450mm(D) x max. 350mm(H)
- Weight: 20kg

Ordering information

- **xrweigh**: Automatic weighing flux
- **3064130**: xrWC Vial & Deep Blue Cap - 250 vials

Shipping dimensions

- Dimensions: 500mm(L) x 500mm(W) x 450mm(H)
- Weight: 50kg

